

CLIP INJECT & HOMEMADE CARTRIDGE



SD Branded product

Single-use, non preloaded intraocular implant injector

TECHNICAL DATA

PARTS & RAW MATERIAL

Body: MABS
Pusher: POM
Sleeve: SEBS
Cartridge: PP
Metal spring

Latex free - DEHP free

Compatible with hydrophilic or hydrophobic lenses

Supplied with viscoglide cartridges:

- 1.8 mm: micron incision
- 2.2 mm: incision size = 2.5 mm

STERILIZATION: Ethylene Oxide

CUSTOMIZATION: color and logo

PACKAGING: 50 devices in blister



The intraocular lens injector enables ophthalmic implant to be safely injected to replace an opacified lens during cataract operations in adults.

Read instructions carefully.

Classe IIa medical device.

CE0459



135 route Neuve
69540 Irigny
France

33 (0) 4 72 39 74 12
irigny@symatesedevise.com
symatese-device.com