

SYRINGE 1ML

SD Branded product



Extra-long 1ML syringe for injecting high-viscosity products such as wrinkle fillers

TECHNICAL DATA

SMOOTH INJECTION

"long design" and specific siliconization

HIGH AESTHETIC & STRENGTH

Glass-like transparency
High resistance to shock and breakage
100% controlled

HIGH CHEMICAL RESISTANCE & EXCELLENT BARRIER PROPERTIES

Prefillable syringes (chemically insert)
Low oxygen permeability
Low water vapor permeability

STERILIZATION: GAMMA

2 options: sterilized / not sterilized

CUSTOMIZATION

Four colors already available
in short time without any investment.
Possibility to develop new colors
throughout short validation and
light investment.

PACKAGING: 100 units/nest in tub



RAW MATERIAL

Injected parts

Syringe barrel: COP

Screw ring: PP

Plunger rod: PC

Finger grip: PC

Components

Tip cap: bromobutyl

Plunger stopper 0.5 ml :

FluroTec WEST réf: 2342

Lubrication: silicone oil



135 route Neuve
69540 Irigny
France

33 (0) 4 72 39 74 12
irigny@symatesedevise.com
symatese-device.com